

```

% Follow along this script with the instructions: "FLIP Algorithm Introduction"

path = '\\millerdata.mit.edu\\amajor\\Earl Miller Lab\\MANUSCRIPTS\\Spectrolaminar pattern\\Minor Revisions_10192023\\2023-11-02 code_AM\\';
file = 'data.mat';

% disp('Select the data.mat file from the directory in which it is saved')
% pause(1);
% [file, path] = uigetfile;
load([path file])

##### Example 1: FLIP Using Default Values
probel = squeeze(data.relpow(863, :, :));
probel(any(isnan(probel), 2), :) = []; % remove NaN rows
disp(['power map is size: ' num2str(size(probel))]);
laminaraxis = 0:0.1:3;
freqaxis = 1:250;
setfreqbool = 1;

% run FLIP with default frequency bands
[startinglowfreq, endinglowfreq, startinghighfreq, endinghighfreq, goodnessvalue, superficialchannel, deepchannel, highfreqmaxchannel, lowfreqmaxchannel, crossoverchannel] = ...
    FLIPAnalysis(probel, laminaraxis, freqaxis, setfreqbool);

##### Example 2: FLIP Using User-Defined Frequency Bin Values
lowfreqrange = [1 30];
highfreqrange = [100 140];
[startinglowfreq, endinglowfreq, startinghighfreq, endinghighfreq, goodnessvalue, superficialchannel, deepchannel, highfreqmaxchannel, lowfreqmaxchannel, crossoverchannel] = ...
    FLIPAnalysis(probel, laminaraxis, freqaxis, setfreqbool, lowfreqrange, highfreqrange);

##### Example 3: Finding Optimal Frequency Bin Values with vFLIP
probe2 = squeeze(data.relpow(807, :, :));
probe2(any(isnan(probe2), 2), :) = [];
setfreqbool = 0;
[startinglowfreq, endinglowfreq, startinghighfreq, endinghighfreq, goodnessvalue, superficialchannel, deepchannel, highfreqmaxchannel, lowfreqmaxchannel, crossoverchannel] = ...
    FLIPAnalysis(probe2, laminaraxis, freqaxis, setfreqbool);

##### Example 4: Generating the spectrolaminar pattern (a.k.a. relative power map) from example LFP data #1
% NOTE: fieldtrip package is required for code beyond this point
% https://www.fieldtriptoolbox.org/download/
% https://www.fieldtriptoolbox.org/faq/should_i_add_fieldtrip_with_all_subdirectories_to_my_matlab_path/
% addpath('/your/directory/for/fieldtrip-20231025');
addpath('\\millerdata.mit.edu\\common\\Alex\\fieldtrip-20191111;');
ft_defaults;
global ft_default
ft_default.showcallinfo = 'no';
ft_default.trackusage = 'no';
ft_warning off;

lfp1 = data.example1_vlPFC_lfp;
[relpow1] = relpow_from_rawLFP(lfp1); % lfp should be in dimension format: nchans x trialltime x ntrials % assumes lfp data in 1 kHz resolution
% requires fieldtrip toolbox
FLIPAnalysis(relpow1, 0:size(relpow1,1)-1, 1:size(relpow1,2), 1); % plot

##### Example 5: Generating the spectrolaminar pattern (a.k.a. relative power map) from example LFP data #2
lfp2 = data.example2_7A_lfp;
[relpow2] = relpow_from_rawLFP(lfp2);
FLIPAnalysis(relpow2, 0:size(relpow2,1)-1, 1:size(relpow2,2), 1); % plot

##### Example 6: Plotting Sample Normalized Power Data from any probe in relpow matrix
row = 752; % there are 942 rows (942 probes in data.relpow) % pick any row
relpow_any_row = squeeze(data.relpow(row, :, :));
relpow_any_row(any(isnan(relpow_any_row), 2), :) = [];
FLIPAnalysis(relpow_any_row, 0:size(relpow_any_row,1)-1, 1:size(relpow_any_row,2), 1); % plot

##### Example 7: Generate Area-average Spectrolaminar Patterns
brain_area = 'MST';
brain_area_num = 2;
index = find([data.meta.brain_area_num]' == brain_area_num);
mean_relpow = squeeze(nanmean(data.relpow(index, :, :), 1));
FLIPAnalysis(mean_relpow, 0:size(mean_relpow,1)-1, 1:size(mean_relpow,2), 1); % plot
sgtitle(['mean LFP relative power (n = ' num2str(length(index)) ' )']);

##### Example 8: Plot Current Source Density for any probe

```

```

row = 752;
CSD = squeeze(data.CSD(row,:,:));
CSD(all(isnan(CSD),2),:) = []; % remove NaN rows
x = figure; imagesc(CSD);set(gca, 'YDir', 'reverse');
ylim([1 size(CSD,1)]); title('CSD'); xlabel('Time (ms)'); ylabel('Channel');
cb=colorbar; ylabel(cb, 'Normalized CSD');

function [relpow] = relpow_from_rawLFP(lfp)

% lfp should be in format nchans x trialtime x ntrials
% assumes lfp data in 1 kHz resolution
% require fieldtrip toolbox

foi = 1:150; % frequencies of interest

nchannels = size(lfp, 1); ntrials = size(lfp, 3); ntrialtime = size(lfp, 2);
% prep ft format
dataLFP = [];
for t = 1:ntrials
    dataLFP.time{t,1} = 0.001:0.001:ntrialtime*0.001; % in seconds
    dataLFP.trial{t,1} = zeros(nchannels, ntrialtime); %nchannels x ntrialtime
    dataLFP.trial{t,1} = lfp(:, :, t);


end
dataLFP.fsample = 1000;
dataLFP.label = {};
%make each channel unique
for c = 1:nchannels
    dataLFP.label{c} = ['ch' num2str(c)];
end
dataLFP.trialinfo = zeros(ntrials, 1);
dataLFP.trialinfo(:, 1) = 1:ntrials;

cfg = [];
cfg.method = 'mtmfft';
cfg.taper = 'hanning';
cfg.output = 'pow';
cfg.keeptrials = 'yes';
cfg.foi = foi;
cfg.pad = 'nextpow2';
pow = ft_freqanalysis(cfg, dataLFP);
meanpow = squeeze(mean(pow.powspectrm));
relpow = meanpow ./ max(meanpow);
end


```

power map is size: 31 250


Default Frequency Bin Values. $G = 1.6792$




Set Frequency Bin Values. $G = 1.4549$




vFLIP Crossover. $G = 1.6711$




Default Frequency Bin Values. $G = 0.78219$



Default Frequency Bin Values. $G = 0.85214$



Default Frequency Bin Values. $G = 1.0753$



mean LFP relative power (n = 145)

